

NEW

ANCHOR
DIAGNOSTICS

ANCHOR BEAD 32

Ultra-Fast Nucleic Acid Purification



The ANCHOR BEAD 32 system delivers fully automated extraction and concentration of pathogen DNA or RNA, powered by SpeedPURE magnetic-bead technology for fast and efficient nucleic acid recovery.


SMARTER EXTRACTION. FASTER ANSWERS.

Universal
protocol




plasma, whole
blood, stool, urine,
CSF, swabs, ...

Flexible
throughput



1–32
samples per run

Compact
design



42 × 41 cm
benchtop
footprint


Efficiency




standardized
workflow,
prefilled cartridges

Simplify your workflow, shorten turnaround times, and deliver actionable results faster. Standardized protocols using SpeedPURE reagents reduce complexity, while flexible throughput ensures your lab can process exactly what's needed - without delay.

SpeedPURE Technology - Ready-to-Use Extraction Reagents



SpeedPURE
Go Cartridge for
single samples




SpeedPURE
Multi Plates for
sample batches

- Pre-filled & sealed SpeedPURE Go Cartridge for individual samples
- Pre-filled & sealed SpeedPURE Multi Plate for batch runs
- Reduced hands-on time and contamination risk

From Sample to Report in 60 Minutes

Unlock peak workflow speed with ANCHOR extraction and PCR kits.

Accelerate molecular pathogen detection with a compact system that extracts PCR-ready DNA and RNA in just 15 minutes. Combined with 30-minute ANCHOR PCR Kits, your lab can deliver results in under 60 minutes - faster reporting with greater efficiency.



< 60 min

Ultra-Fast Turnaround

- ~15 min extraction on ANCHOR BEAD 32 platform
- 30 min PCR with ANCHOR PCR Kits
- Sample-to-result < 60 min
- Improve flexibility in lab operations
- High-speed workflow ensures timely results

SIMPLIFIED. STANDARDIZED. ACCELERATED.

The ANCHOR BEAD 32 instrument is intended for Research Use Only (RUO). Not for use in diagnostic procedures.

Contact us today to request your demo or learn more about transforming your lab's workflow.



ANCHOR Diagnostics GmbH
Grandweg 64
22529 Hamburg, Germany
phone: +49 40 52 471 62 0
www.anchor-diagnostics.com