

WHY THIS MATTERS FOR YOUR LABORATORY:

- Report your PCR results significantly faster
- Process any sample without delay
- Increase throughput and flexibility
- Enhance staff efficiency and eliminate delays
- Maximize utilization of existing PCR cyclers
- Add speed and flexibility to your existing molecular platforms

Made for
Speed & Flexibility



Because Speed Matters.



ANCHOR Diagnostics GmbH
Grandweg 64, 22529 Hamburg, Germany
phone: +49 40 52 471 62 0
email: info@anchor-diagnostics.com
anchor-diagnostics.com

Precision Molecular Pathogen Detection – Results in Under 60 Minutes

Simple. Standardized. Ultra-Fast.

HIGH-SPEED FOR YOUR MOLECULAR WORKFLOW



1 Any Sample. One Standardized Workflow.

Harmonized preparation.
Consistent performance.

- One standardized workflow for diverse sample materials
- Whole blood, plasma, CSF, swabs, urine, stool
- Consistent performance across matrices

2 SpeedPURE® – Reagents for Speed DNA and RNA Extraction.*

Eliminate preparation steps.
Start immediately.

- Optimized for ~15 min extraction
- Pre-filled and sealed for instant use
- Universal DNA & RNA protocol for diverse sample types
- Three kit formats for every throughput need

3 ANCHOR BEAD 32 – Ultra-Fast Nucleic Acid Extraction*

Engineered for High-Speed Routine Extraction.

- Fully automated DNA and RNA extraction in ~15 min
- Fully flexible for 1-32 samples
- Preinstalled protocols for different sample volumes
- Integrated UV and HEPA protection

4 ANCHOR PCR Kits – Designed for Speed & Reliability.

Fast results.
Harmonized performance.

- Results in 30 minutes
- Harmonized pipetting and PCR profiles
- Broad PCR cycler compatibility
- Excellent analytical performance

From sample to result in one seamless workflow.
Scan to see how it works:



*SpeedPURE® and ANCHOR BEAD 32 currently available as research use only (RUO). Available as CE IVD version soon. SpeedPURE® is a registered trademark of ANCHOR Diagnostics. All rights reserved.